## RESTYLANE® REFYNE™ INSTRUCTIONS FOR USE

I. COMPOSITION	
Cross-linked hyaluronic acid	20 mg/mL
Lidocaine hydrochloride	3 mg/mL
Phosphate buffered saline pH 7	qs ad 1 mL

### II. DESCRIPTION

Restylane Refyne is a sterile, biodegradable, transparent gel of non-animal cross-linked hyaluronic acid with the addition of lidocaine hydrochloride 3 mg/mL. The gel is supplied in a pre-filled plastic syringe. The contents of the syringe are sterilized using moist heat. The syringe is packaged individually in a blister, with two 30G x 1/2" Ultra thin wall needles. The needles have been sterilized using irradiation. The product is for single use only. To ensure traceability the package includes patient record labels that should be attached to patient records.

### III. INTENDED USE

Restylane Refyne is intended for injection into the mid to deep dermis for correction of moderate facial wrinkles and folds such as nasolabial folds. Lidocaine is added to the formulation to diminish the pain resulting from the injection during the treatment.

The product is intended to be used only by authorized personnel in accordance with local legislation, trained in the appropriate injection techniques. Before the first treatment session, it is recommended to contact your local Galderma representative or Restylane distributor for information about training opportunities.

### IV. MODE OF ACTION

The product adds volume to the tissue thereby restoring the skin contours of the face to the desired level. The volume and lifting capacity originate from the ability of cross-linked hyaluronic acid to bind water.

### V. CONTRAINDICATIONS

- Do not use in patients with a history of hypersensitivity to streptococcal proteins, as the product may contain trace amounts of such material.
- Do not use in patients with severe allergies manifested by a history of anaphylaxis or history or presence of multiple severe allergies.
- Do not use in patients with known hypersensitivity to lidocaine or to amide-type local anaesthetics.
- Patients presenting with porphyria

### VI. WARNINGS

- Use at specific sites where there is active disease, such as inflammation (skin eruption such as cysts, pimples, rashes or hives), infection or tumours, in or near the intended treatment site should be avoided until the underlying process has been controlled.
- This product must not be injected intramuscularly or intravascularly.
- Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.
- Localized superficial necrosis and scarring may occur after injection in or near vessels, such as the nose and glabella area. It is thought to result from the injury, obstruction, or compromise of blood vessels. Special caution should be taken if the patient has undergone a prior surgical procedure in the planned treatment area. Areas with limited collateral blood flow has an increased risk of ischemia. Aspiration prior to injection is recommended.

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- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.
- Patients with bleeding disorders or patients using substances that affect platelet function, thrombolytics or anticoagulants may as with any injection, experience increased bruising or bleeding at the injection site.
- This product should not be mixed with other products prior to injection.

### VII. PRECAUTIONS

- In order to minimize the risks of potential complications (perforation or compression of vessels, nerves and other vulnerable structures), this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.
- Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.
- Patients with unattainable expectations are not suitable candidates for treatment.
- The safety for use during pregnancy, in breastfeeding females or in patients under 18 years has not been established. Do not use in children.
- Injection procedures are associated with a risk of infection. Aseptic technique and standard practice to prevent cross-infections are to be followed.
- Do not use the product if the package is open or damaged, or if the expiry date or lot number is illegible.
- Do not use the product if the contents of the syringe are cloudy.
- This product should be used with caution in patients on immunosuppressive therapy.
- This product should be used with caution in patients with a tendency to form hypertrophic scars or any other healing disorders.
- At each treatment session, a maximum dosage of 2 mL per treatment site is recommended.
- Avoid injecting into areas in close proximity to permanent implants, as this could aggravate latent adverse events or interfere with the aesthetic outcome of the treatment. Limited data is available on injecting into an area where a non-permanent implant other than hyaluronic acid has been placed.
- Injection too superficially, or in facial areas with limited soft tissue support or soft tissue cover, or thin skin, such as the periorbital area, may result in contour irregularities and palpable lumps and/or bluish discoloration.
- Patients with pre-existing pigmented dark lower eye lid circles, thin skin and pre-existing tendency toward oedema formation are not suitable candidates for treatment of the lower periorbital region as it may be associated with prominent discoloration and excessive swelling due to fluid build-up.
- Patients should avoid excessive sun, UV lamp exposure and extreme temperatures until any initial swelling and redness has resolved.

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If laser treatment, chemical peeling or any other procedure based on active dermal response is performed after treatment with this product there is a theoretical risk of eliciting an inflammatory reaction at the implant site. This also applies if the product is administered before the skin has healed completely after such a procedure.

- Injection procedures can lead to reactivation of latent or subclinical herpes viral infections.
- Post-inflammatory pigmentation changes may occur after dermal filler injections in people with dark skin (Fitzpatrick Type IV-VI).
- Individual variation and treatment area may affect the bio-degradation of this product, in rare cases product remnants have been detected in tissue when the clinical effect has returned to baseline.
- This product is packaged for single use. Do not re-sterilize.
- Considerations should be given to the total dose of lidocaine administered if dental block or topical administration of lidocaine is used concurrently. High doses of lidocaine (more than 400 mg) can cause acute toxic reactions manifesting as symptoms affecting the central nervous system and cardiac conduction.
- Lidocaine should be used with caution in patients receiving other local anaesthetics or agents structurally related to amide-type local anaesthetics e.g., certain anti-arrhythmics, since the systemic toxic effects can be additive.
- Lidocaine should be used cautiously in patients with epilepsy, impaired cardiac conduction, severely impaired hepatic function or severe renal dysfunction.

### VIII. ADVERSE EVENTS

Patients must be informed of the potential risks and adverse events related to the injection procedure and to the use of this product.

### Anticipated injection-related reactions

Injection-related reactions are expected to occur after treatment with Restylane Refyne. These include but are not limited to injection-related reactions such as bruising, erythema, itching, swelling, pain and tenderness at implant site. Typically, resolution is spontaneous within one week after injection.

### Post marketing adverse event reporting

The following post marketing adverse events have been reported from worldwide sources after treatment with Restylane Refyne (non-exhaustive list). The most commonly reported included transient swelling/edema and with immediate onset and onset up to several weeks after treatment.

The following events were also reported: Paps/induration, Erythema, Bruising/bleeding, Pain/tenderness, Papules/nodules, Device ineffective, Discoloration, Infection/ abscess including pustule and purulent discharge, Ischemia/necrosis including pallor, ulcer and vascular occlusion, Inflammation, Hypersensitivity/angioedema, Other injection site reactions and skin reactions including burning sensation, exfoliation, dryness and warmth, Device dislocation, Pruritus, Eye disorders including eye swelling, ocular discomfort and visual impairment such as blurred vision, Deformity/asymmetry including overcorrection, Blisters/vesicles, Neurological symptoms including hypoesthesia and paraesthesia, Granuloma/foreign body reaction, Scar/skin atrophy, Rash, Discharge, Acne, Capillary disorders such as telangiectasia, Reactivation of herpes infection, Urticaria, Dermatitis, Muscle disorder including muscle twitching. Non-dermatological events including chills, dizziness, dyspnea and headache, Other dermatological events such as chapped lips.

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Vascular compromise may occur due to an inadvertent intravascular injection or as a result of vascular compression associated with implantation of any injectable product. This may manifest as blanching, discoloration, necrosis or ulceration at the implant site or in the area supplied by the blood vessels affected, or rarely as ischemic events in other organs due to embolization. Rare but serious cases of ischemic events associated with temporary or permanent vision impairment, blindness, cerebral ischemia or stroke have been reported following facial aesthetic treatments. Symptoms of inflammation at the implant site commencing either shortly after injection or after a delay of up to several weeks have been reported. In case of unexplained inflammatory reactions infections should be excluded and treated if necessary since inadequately treated infections may progress into complications such as abscess formation. Treatment using only oral corticosteroids without concurrent antibiotic treatment is not recommended. The prolonged use of any medication, e.g., corticosteroids or antibiotics in treatment of adverse events has to be carefully assessed, since this may carry a risk for the patient. In case of persistent or recurrent inflammatory symptoms, consider removal of the product by aspiration/drainage, extrusion or enzymatic degradation (use of hyaluronidase has been described in scientific publications). Before any removal procedure is performed, the swelling may be reduced by using NSAID for 2-7 days or a short course of corticosteroids for less than 7 days, in order to more easily palpate any remaining product. For patients who have experienced clinically significant reactions, a decision for retreatment should take into consideration the cause and significance of previous reactions. For reporting of adverse events contact your local Galderma representative or distributor for this product.

### IX. PERFORMANCE

In a randomized controlled multicenter study with Restylane Refyne for correction of moderate to severe nasolabial folds, 78.0% and 62.2% of subjects maintained a clinically significant improvement (>1 grade improvement in Wrinkle Severity Rating Scale [WSRS] score compared to baseline) at 36 and 48 weeks after treatment, respectively.

### X. ASSEMBLY OF NEEDLE TO SYRINGE

It is important to use a sterile, appropriate needle or blunt cannula. Suitable needles (30G x 1/2" Ultra thin wall) are supplied with the syringe in the blister pack. As an alternative, a blunt thin walled cannula with a recommended size of 27G can be used. Use surgical gloves, remove the cap from the needle and unscrew the tip cap from the syringe. Hold firmly around the syringe barrel and grasp the needle shield with the other hand. Screw the needle tight onto the syringe by simultaneously pushing and rotating firmly until the needle is completely locked. To ensure proper assembly, minimize the gap between the needle shield and the syringe. See Figure 1. Remove the needle shield just before injection by pulling it straight out. Do not rotate. Note! Improper assembly may cause leakage or needle disconnection. Assembly instruction applies to both needle and blunt cannula.

### XI. TREATMENT PROCEDURE

- Inform the patient about the treatment procedure and expected results. Ensure the patient receives appropriate post-treatment care instructions and is made aware of signs and symptoms of potential complications.
- Cleanse the area to be treated with an antiseptic and allow it to dry before injection.

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- To avoid breakage of the needle/cannula, do not attempt to bend it before or during treatment. If the needle gets bent, discard it and complete the procedure with a replacement needle.
- Before injecting the product, depress the plunger rod carefully until a small droplet is visible at the tip of the needle/cannula.
- Align the bevel of the needle by turning the syringe on its axis.
- If a blunt cannula is used, an entry point is made in the skin, for example with a sharp needle of appropriate size. During injection, keep the side hole of the cannula facing downwards away from the skin surface.
- Aspiration is recommended prior to injection in order to reduce the risk of inadvertent injection into a blood vessel.
- Inject the gel slowly by gently pressing down on the plunger rod with the thumb or palm of the hand. Do not apply excessive pressure to the syringe at any time. Presence of scar tissue may impede advancement of the needle/cannula. If resistance is encountered the needle/cannula should be partially withdrawn and repositioned or fully withdrawn and checked for function.
- Choose from a variety of injection techniques, i.e. serial puncture, linear threading or cross-hatching. It is recommended to change needle/cannula for each new treatment site.
- Defects should be fully corrected, but not overcorrected, at each treatment session. Gently massage the treated area after injection.
- If "blanching" of the skin is observed the injection should be stopped immediately and the whitened area should be massaged gently until it returns to a normal colour before continuing with the injection.
- If the treated area is swollen directly after the injection, an ice pack with adequate protective cloth can be applied on the site for a short period. Ice should be used with caution if the area is still numb from anaesthetic to avoid thermal injury.
- Additional treatments may be necessary to achieve and maintain the desired level of correction.

The syringe, disposable needle/cannula and any unused material must be discarded immediately after the treatment session and must not be reused due to the risk for contamination of the unused material and the associated risks including infections. Dispose should be in accordance with accepted medical practice and applicable national, local or institutional guidelines.

### XII. SHELF LIFE AND STORAGE

Do not use after the expiry date indicated on package. Store up to 25°C. Protect from freezing and sunlight.

### XIII. MANUFACTURER

Q-Med AB, Seminariegatan 21, SE-752 28 Uppsala, Sweden  
Phone +46(0)18 474 90 00, Fax +46(0)18 474 90 01  
www.galderma.com, e-mail: info.q-med@galderma.com

### XIV. DISTRIBUTOR

Galderma Canada Inc.  
Thorhill, ON L3T 7V9  
1 800 467-2081

Restylane and Galderma are registered trademarks.

Manufacturer of the needles in the blister pack:

TSK LABORATORY, JAPAN,  
2-1-5, Hiraynagi-cho, Tochigi-shi,  
Tochigi-ken, 328-0012 - JAPAN

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Injectable Gel with Lidocaine

## Assembly instruction/Instructions d'assemblage

For further information, see section X.

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