## Restylane® Silk – Instructions for Use

**Composition**  
Hyaluronic acid, stabilized 20 mg/mL  
Lidocaine hydrochloride 3 mg/mL  
Phosphate buffered saline pH 7

**Description**  
Restylane® Silk is a sterile, transparent hyaluronic acid gel supplied in a glass syringe together with a 30G or 29G Thin Wall (TW) needle. This is the single use only syringe and is a unique formulation, stabilized hyaluronic acid (NASHA™). Hyaluronic acid is a natural polysaccharide which occurs as an important structural element in the extracellular matrix of connective tissues as well as in the synovial fluid and bone. Hyaluronic acid belongs to a group of very few substances which are identical in all living organisms.

**Mode d'action**  
Restylane® Silk acts by adding volume to the tissue, thereby restoring skin contours to the desired level of correction. Restylane® Silk is a temporary filler that will remain in tissue under isovolumic degradation.

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Holding the syringe by its ribbed part (C) of the white closure system (base/stopper).

With your other hand take hold of the white cap (A) at the end of the closure system and gently pull and forth carefully until disengaged and can be pulled off (cap A will be broken).

Do not reuse.

Do not touch the syringes tip (B) to keep it sterile.

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Assemble the needle to syringes

Use the thumb and forefinger to hold firmly around both the glass syringe barrel and the luer-lock adapter. Grasp the needle shaft with the other hand. To facilitate proper assembly, both push and rotate firmly.

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Assembly of needle to syringes

Use the thumb and forefinger to hold firmly around both the glass syringe barrel and the luer-lock adapter.

Grasp the needle shaft with the other hand. To facilitate proper assembly, both push and rotate firmly.

**Contraindications**  
Patients known to hypersensitivity to hyaluronic acid, lidocaine or amide type local anaesthetics. Patients with bleeding disorders.

**Special caution**  
• Special caution is required to avoid perforation or compression of vessels, nerves and other vulnerable structures.

• Localized scleroderma and scleritis may occur with lidocaine. Special caution should be taken if the patient has a history of scleroderma or scleritis. Risk of scleroderma or scleritis is increased if the patient has a history of Eosinophilic granuloma, sarcoidosis, vasculitis, cerebral ischemia or cerebral hemorrhage, leading to stroke, sinus or nasal damage, and/or peripheral nerve palsies.

• Rare but serious adverse associated with the intravascular injection of soft tissue fillers in the face have been reported and include transient facial palsy, transient facial nerve palsy, trigeminal neuralgia, cerebral ischemia or cerebral hemorrhage, leading to stroke, sinus or nasal damage, and/or peripheral nerve palsies.

• Patients with preexisting pigmentary dark skin, including Fitzpatrick skin types IV-VI, are at increased risk for transient or permanent hypoesthesia or anesthesia after injection.

• Patients with preexisting skin conditions or diseases associated with skin changes are not suitable candidates for treatment.

• Do not use if the product is damaged.

• Specific considerations relevant to the use of this product:

• Do not mix with other products prior to injection or device.

• Rare but serious adverse products such as inflammation, infection or tumours can occur at the injection site.

• Patients should avoid excessive sun or extreme heat, as well as any initial swelling or redness due to the injection.

• If laser treatment, chemical peeling or any other procedure causing significant tissue damage is performed, wait at least 2 weeks before injecting Restylane® Silk.

• Health care practitioners are encouraged to minimize the risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

• Hypersensitivity reactions are associated with a risk of infection. Aspirate technique and standard practice to prevent cross-infections must be observed.

• Lidocaine should be used with caution in patients receiving other local anaesthetics or agents structurally related to amide-type lidocaine, e.g. cloba-antiarrhythmics, since the systemic toxic effects can be additive.

**Precautions**  
General considerations relevant to medical devices

• To minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience and who have knowledge of the indications, contraindications and side effects of the product.

• Health care practitioners are encouraged to minimize the risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

• Considerations should be given to the total dose of lidocaine administered if dental block or topical anaesthesia of lidocaine is used concomitantly with Restylane® Silk. Lidocaine (up to 400mg) can cause acute toxic reactions as manifested by symptoms affecting the central nervous system.

• Lidocaine should be used with caution in patients receiving other local anaesthetics or agents structurally related to amide-type lidocaine, e.g. cloba-antiarrhythmics, since the systemic toxic effects can be additive.

Other potential adverse events that have been reported following injection of other Restylane products and that may occur with using Restylane® Silk include: localised pain, swelling, bruising, haematomas, Purpura, papule/nodules, Inflammation, Eosinophilia, Skin rash, loss of colour, effect of injection and/or bruising formation, Drowsiness, Hyperactivity, Neurogical symptoms, including paraesthesia, Eye pain, Granuloma.

Other adverse events that have been reported following injection of other Restylane products and that may occur with using Restylane® Silk include: localised pain, swelling, bruising, haematomas, Purpura, papule/nodules, Inflammation, Eosinophilia, Skin rash, loss of colour, effect of injection and/or bruising formation, Drowsiness, Hyperactivity, Neurogical symptoms, including paraesthesia, Eye pain, Granuloma.

Other adverse events that have been reported following injection of other Restylane products and that may occur with using Restylane® Silk include: localised pain, swelling, bruising, haematomas, Purpura, papule/nodules, Inflammation, Eosinophilia, Skin rash, loss of colour, effect of injection and/or bruising formation, Drowsiness, Hyperactivity, Neurogical symptoms, including paraesthesia, Eye pain, Granuloma.

Other adverse events that have been reported following injection of other Restylane products and that may occur with using Restylane® Silk include: localised pain, swelling, bruising, haematomas, Purpura, papule/nodules, Inflammation, Eosinophilia, Skin rash, loss of colour, effect of injection and/or bruising formation, Drowsiness, Hyperactivity, Neurogical symptoms, including paraesthesia, Eye pain, Granuloma.

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